

Manual visual inspection systems.

- ▣ Accurate and rapid inspection of all container and product types
- ▣ Meet all global standards and regulations for visual inspection systems, including those of the World Health Organization, United States, European and Japanese pharmacopeias + CE certificated
- ▣ Non-glare black and white background + Two (2) fluorescent lamps or optional LED lighting
- ▣ Utilities: 230 V, 1 phase, 50 Hz

MIH-PORT

- ▣ Dimensions: 708 mm (W) x 551 mm (H) x 505 mm (D)
- ▣ Black powder coat finish
- ▣ White base with padded elbow rests

Advantages of MIH-Port

- + Compact handling and storage - foldable and portable
- + Elbow rests for operator's comfort



MIH-DX

- ▣ Dimensions: 686 mm (W) x 616 mm (H) x 544 mm (D)
- ▣ Stainless steel enclosure
- ▣ Adjustable position angle
- ▣ Courtesy side shield

Advantages of MIH-DX

- + Ideal light conditions due to side shields
- + Available with tyndall light to easier detect particles in liquids



MIH-LX

- Dimensions: 737 mm (W) x 519 mm (H) x 470 mm (D)
- Black powder coat finish
- Sensor verification of hand (and container) in the inspection zone with programmable inspection/pacing time
- White base with padded elbow rests

Advantages of MIH-LX

- + Height adjustable lighting and background
- + Equipped with timer, two LED indicator lamps, two position sensitive infra-red photoelectric sensors, one alarm buzzer, one "Timer Reset" button and an operator interface panel



MIH-1

- Dimensions: 1524 mm (W) x 1527 mm (H) x 692 mm (D)
- Extruded aluminum frame and black plastic panels
- White Corian® work surface
- Adjustable side trays or tube holders (syringes)
- Ergonomic lift system with adjustable height range of 30 to 46 in. (76 to 117 cm) for standing or sitting inspection

Advantages of MIH-1

- + Electrically operated lift mechanism to raise hood 16" to facilitate ergonomic sitting and standing positions
- + Side shields for better light conditions



Additional options	MIH-Port	MIH-DX	MIH-LX	MIH-1
Fluorescent lighting	yes	yes	yes	yes
Fixed white LED lighting	yes	yes	yes	yes
Adjustable intensity white LED lighting	yes	yes	yes	yes
Magnification lens - 1.75X		yes	yes	yes
Tyndall lighting		yes		yes
Automatic timing feature			yes	yes
Vial tray holders				yes
Syringe tub holders				yes

Connect with our team!

Syntegon Technology GmbH

packaging-ph@syntegon.com

www.syntegon.com