



**Manual visual
inspection systems.**

Manual visual

inspection systems.

- Accurate and rapid inspection of all container and product types
- Meet all global standards and regulations for visual inspection systems, including those of the World Health Organization, United States, European and Japanese pharmacopeias + CE certificated
- Non -glare black and white background with various LED overhead lighting options
- Utilities: 100 - 230 VAC, 50/60 Hz, 1 phase

MIH-PORT



- Dimensions: 708 mm (W) x 551 mm (H) x 505 mm (D)
- Black powder coat finish
- White base with padded elbow rests

Advantages of MIH-Port

- + Compact handling and storage - foldable and portable
- + Elbow rests for operator's comfort

MIH-DX




- Dimensions: 686 mm (W) x 616 mm (H) x 544 mm (D)
- Stainless steel enclosure
- Adjustable angled background
- Courtesy side shield

Advantages of MIH-DX

- + Ideal light conditions due to side shields
- + Available with tyndall light to easier detect particles in liquids

MIH-DL



- ❑ Dimensions in mm: 692 (W) x 538 (H) x 432 (D)
- ❑ Matte black aluminum construction minimizes glare
- ❑ Dual adjustable LED lighting
- ❑ Courtesy side shield

Advantages of MIH-DL

- + Optimal environment for difficult to inspect (DIP) products
- + LED lighting features isotropic lenses to ensure homogenous lighting with minimal specular reflection or hotspots

MIH-LX



- ❑ Dimensions: 737 mm (W) x 519 mm (H) x 470 mm (D)
- ❑ Black powder coat finish
- ❑ Sensor verification of hand (and container) in the inspection zone with programmable inspection/pacing time
- ❑ White base with padded elbow rests

Advantages of MIH-LX

- + Height adjustable lighting and background
- + Equipped with timer, two LED indicator lamps, two position sensitive infra-red photoelectric sensors, one alarm buzzer, one "Timer Reset" button and an operator interface panel

MIH-IV



- ❑ Dimensions in mm: 768 (W) x 703 (H) x 491 (D)
- ❑ Optimal inspection environment for IV bags
- ❑ Large 457 x 508 mm background
- ❑ Accommodates up to 2 liter IV bags

Advantages of MIH-IV

- + A sliding hook enables the operator to hang the IV bag for clear direct inspection
- + Sliding the bag promotes the agitation of particulate movement enhancing inspection capability

MIH-1



- Dimensions: 1524 mm (W) x 1527 mm (H) x 692 mm (D)
- Extruded aluminum frame and black plastic panels
- White Corian® work surface
- Adjustable side trays or tub holders (syringes)
- Ergonomic lift system with adjustable height range of 30 to 46 in. (76 to 117 cm) for standing or sitting inspection

Advantages of MIH-1

- + Electrically operated lift mechanism to raise hood 16" to facilitate ergonomic sitting and standing positions
- + Side shields for better light conditions

MIH-1S



- Dimensions in mm: 914 (W) x 1524 (H) x 686 (D)
- Ergonomic lift system with adjustable countertop, height range of 76 to 97 cm

Advantages of MIH-1S

- + Electrically operated lift mechanism to facilitate ergonomic working position
- + Side shields for better light conditions

Lighting and machine options.

ADDITIONAL OPTIONS	MIH-Port	MIH-DX	MIH-DL	MIH-LX	MIH-1	MIH-1S
Single light intensity LED lighting	yes	yes		yes	yes	yes
Dual light intensity LED lighting	yes	yes		yes	yes	yes
Adjustable intensity white LED lighting	yes	yes	yes	yes	yes	yes
Tyndall lighting		yes			yes	yes
Magnification lens - 1.75X		yes	yes	yes	yes	yes
Automatic timing feature				yes	yes	yes
Vial tray holders					yes	
Syringe tub holders					yes	

Connect with our team!

Syntegon Technology GmbH

packaging-ph@syntegon.com | www.syntegon.com